

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

Page 1

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|-------------------------------|-------------------------------------|---|--------------------------------|-------------|
| Ecolab® | A-456-II | | 0.4% A-256-II (1:256 Water) | 02-Dec-14 |
| Ecolab® | Actichlor Plus | Sodium Dichloroisochyanurate 30-60% | 0.2% Actichlor Plus (2000PPM) | 20-Nov-12 |
| Quantum Pacific Limited | Actimax | | According to Manufacturer Spec | 11-Apr-06 |
| SCIENCECompany | Ammonium Chloride | Cocoalkaydimethylbenzyl ammonium chloride | According to Manufacturer Spec | 17-May-19 |
| Milton® | Anti-Bacterial Tablets | Sodium Dichloroisochyanurate 30-60% | 1 Tablet:2L Water | 19-Jul-12 |
| Ecolab® | Asepti Active Liquid | | According to Manufacturer Spec | 20-Nov-12 |
| Ecolab® | Asepti Active Wipes | Dual-Quaternary/Alcohol | Ready to Use | 20-Nov-12 |
| Synergy Health | Azomax™ Multi Surface Wipes | Benzylnoniumchloride 21.6% | Ready to Use | 15-Mar-13 |
| Synergy Health | Azowipe™ Hard Surface Wipes | Benzylnoniumchloride 2% | Ready to Use | 15-Mar-13 |
| PreVasive® | BAC Botanical Antimicrobial Cleaner | Thyme Oil | According to Manufacturer Spec | 22-May-19 |
| Spartan Chemical Company Inc. | BNC-15™ | | According to Manufacturer Spec | 11-Aug-16 |
| Metrex® Research | CaviCide® | Isopropanol | According to Manufacturer Spec | 06-Oct-16 |
| Metrex® Research | CaviCide1™ | Isopropanol | According to Manufacturer Spec | 05-Jan-17 |
| Metrex® Research | CaviWipes1™ | Isopropanol | Ready to Use | 05-Jan-17 |
| Metrex® Research | CaviWipes™ | Isopropanol | Ready to Use | 06-Oct-16 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

Page 2

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|---------------------------------|---|--|-----------------------------------|-------------|
| Agar Cleaning Systems Pty. Ltd. | Chloradet | Sodium Hypochlorite <10% | 0.1% Cloradet (1:1000 Water) | 06-Jun-16 |
| Apollo Corporation | CID A-L™ II Quaternary Disinfectant | Quaternary Ammonium Compounds | According to Manufacturer Spec | 26-Sep-12 |
| Wexford Labs, Inc. | CleanCide® Wipes | Citric Acid | Ready to Use | 15-May-19 |
| Clorox Disinfecting Wipes | Clorox Disinfecting Wipes | | According to Manufacturer Spec | 14-Feb-11 |
| Alco Mandala Products | DC-7 | | According to Manufacturer Spec | 16-Feb-12 |
| Clinell® | Detergent Wipes | Propane 1,2-diol | According to Manufacturer Spec | 27-May-09 |
| ArjoHuntleigh | Disinfectant Cleanser IV | | According to Manufacturer Spec | 23-Nov-09 |
| Clorox® | Dispatch® Disinfectant | | According to Manufacturer Spec | 16-Feb-12 |
| Tristel Solutions LTD | Distel Surface Disinfectant (Formerly Trigene®) | Polymeric Biguanide Hydrochloride <10% | According to Manufacturer Spec | 24-Oct-06 |
| Sealed Air Diversey Care | Divercleanse | Sodium Hypochlorite <5%, Sodium Peroxide <2% | 3.1% Divercleanse (3.1:100 Water) | 15-Dec-09 |
| Buckeye® International, Inc. | Eco E23 S23 | | 256 to 1 dilution with water | 05-Jan-17 |
| Velocity Chemicals | Eco-Pine | | 5% Ecopine (1:20 Water) | 11-Apr-06 |
| Clorox HEALTHCARE® | Fuzion™ | | According to Manufacturer Spec | 05-Jan-17 |
| Spartan Chemical Company Inc. | GS Neutral Disinfectant #3502 | | According to Manufacturer Spec | 02-Dec-14 |
| Clinell® | Hand & Surface Sanitizer | Benzalkonium Chloride | According to Manufacturer Spec | 27-May-09 |

CARR NEW ZEALAND LIMITED

 Tel: 64 9 274 0869
 6 Ra Ora Drive, East Tamaki, Auckland 2013, New Zealand
 PO Box 204110, Highbrook, Auckland 2161, New Zealand
 info@carrgroup.co.nz

CARR AUSTRALIA PTY LIMITED

 Tel: 61 2 9636 9525
 1 Spireton Place, Pendle Hill, NSW 2145, Australia
 PO Box 335, Pendle Hill, NSW 2145, Australia
 info@carrgroup.com.au

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

Page 3

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|-------------------------------------|--|---|--------------------------------|-------------|
| 3M | HB Quat Disinfectant | Alkyl dimethyl benzyl ammonium chloride | Ready to Use | 06-Sep-19 |
| Spartan Chemical Company Inc. | hdq®C2 | Ammonium Chloride | According to Manufacturer Spec | 30-Jun-17 |
| Wipesplus | Hospital/Surgical Wipes | Ammonium Chloride | Ready to Use | 15-May-19 |
| Clorox HEALTHCARE® | Hydrogen Peroxide Cleaner Disinfectant Wipes | Hydrogen Peroxide | Ready to Use | 12-Apr-19 |
| Apollo Corporation | Hygena™ | Ethylalcohol 1.4% | According to Manufacturer Spec | 25-Sep-12 |
| Kimberly-Clark® | Isowipes® | Isopropanol 30-60% | Ready to Use | 19-Jul-12 |
| Legge Systems | Legcide HB Plus | | According to Manufacturer Spec | ***** |
| Whiteley Medical | Lemex | | According to Manufacturer Spec | 15-Mar-11 |
| Livingstone International (PTY) LTD | LIV-Wipes Jumbo | 70% Isopropyl Alcohol | Ready to Use | 15-Mar-11 |
| Sulco LTD | Mediwipes | | Ready to Use | 27-May-09 |
| Medline Industries Inc. | Micro-Kill™ Q3 Concentrated Disinfectant, Cleaner & Deodorizer | | According to Manufacturer Spec | 01-Feb-18 |
| Medline Industries Inc. | MICRO-KILL™ AF2 Wipes | Ammonium Chloride | Ready to Use | 15-May-19 |
| Medline Industries Inc. | MICRO-KILL™ Bleach-Free, Alcohol-Free Wipes | Ammonium Chloride | Ready to Use | 15-May-19 |
| Medline Industries Inc. | MICRO-KILL™ One Wipes | Ethyl Alcohol | Ready to Use | 15-May-19 |
| Pro Restore® Products | Microban X-590 Institutional Spray | Diisobutylphenoxyethoxy ethyl dimethyl benzyl ammonium chloride | According to Manufacturer Spec | 11-Aug-16 |

CARR NEW ZEALAND LIMITED

Tel: 64 9 274 0869
 6 Ra Ora Drive, East Tamaki, Auckland 2013, New Zealand
 PO Box 204110, Highbrook, Auckland 2161, New Zealand
 info@carrgroup.co.nz

CARR AUSTRALIA PTY LIMITED

Tel: 61 2 9636 9525
 1 Spireton Place, Pendle Hill, NSW 2145, Australia
 PO Box 335, Pendle Hill, NSW 2145, Australia
 info@carrgroup.com.au

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

Page 4

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|--------------------------------|---|-----------------------|--------------------------------|-------------|
| Anti-Slip Anywhere | Multi Kleen | | According to Manufacturer Spec | 15-Mar-11 |
| Certo® | Non Acid Bathroom Cleaner | | According to Manufacturer Spec | 07-May-12 |
| PreVasive® | Noroxycdiff™ | Hydrogen Peroxide | According to Manufacturer Spec | 22-May-19 |
| Ecolab® | Oasis Pro 12 Neutral Cleaner | Alcohols <10% | 5% Oasis Pro 12 (1:20 Water) | 06-Mar-12 |
| Whiteley Medical | Odacon Spray | | According to Manufacturer Spec | 14-Jun-11 |
| OdorXit Naturals | Odor Eliminator Spray | | Ready to Use | 22-Feb-18 |
| Pal® International | PAL® Detergent Wipes | | Ready to Use | 27-May-09 |
| Pal® International | PAL® Disinfectant Wipes | | Ready to Use | 27-May-09 |
| PCS Process Cleaning Solutions | PCS 7000 | Sodium Hypochlorite | According to Manufacturer Spec | 06-Oct-17 |
| PCS Process Cleaning Solutions | PCS Spray & Wipe Hand and Surface Cleaner | | Ready to Use | 06-Oct-17 |
| Pure Green, LLC | Puregreen 24 (Pure Green) | | According to Manufacturer Spec | 30-Jan-09 |
| Reynard Health Supplies | Reynard Neutral Detergent Wipes | | According to Manufacturer Spec | 08-Feb-07 |
| Reynard Health Supplies | Reynard Surface Sanitizing Wipes | | According to Manufacturer Spec | 27-May-11 |
| Sealed Air Diversey Care | Room Care R2L | | According to Manufacturer Spec | 24-Oct-12 |
| STERI-7® | S-7 Xtra Wipes | | Ready to Use | 28-Nov-16 |

CARR NEW ZEALAND LIMITED

 Tel: 64 9 274 0869
 6 Ra Ora Drive, East Tamaki, Auckland 2013, New Zealand
 PO Box 204110, Highbrook, Auckland 2161, New Zealand
 info@carrgroup.co.nz

CARR AUSTRALIA PTY LIMITED

 Tel: 61 2 9636 9525
 1 Spireton Place, Pendle Hill, NSW 2145, Australia
 PO Box 335, Pendle Hill, NSW 2145, Australia
 info@carrgroup.com.au

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

Page 5

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|-------------------------|--|-----------------------|--------------------------------|-------------|
| Alco Mandala Products | Sanni-Rinse | | According to Manufacturer Spec | 16-Feb-12 |
| Clinell® | Sporicidal Wipes | | According to Manufacturer Spec | 27-May-09 |
| Ecolab® | Spring Clean | | 5% Spring Clean (1:20 Water) | 11-Apr-06 |
| Noble Pine Products Co. | Steri-Fab® | Isopropyl Alcohol | According to Manufacturer Spec | 11-Aug-16 |
| sBioMed, LLC | Steriplex SD Disinfectant Sporicide Wipes (Without Bleach) | Hydrogen Peroxide | According to Manufacturer Spec | 20-May-15 |

CARR NEW ZEALAND LIMITED

 Tel: 64 9 274 0869
 6 Ra Ora Drive, East Tamaki, Auckland 2013, New Zealand
 PO Box 204110, Highbrook, Auckland 2161, New Zealand
 info@carrgroup.co.nz

CARR AUSTRALIA PTY LIMITED

 Tel: 61 2 9636 9525
 1 Spireton Place, Pendle Hill, NSW 2145, Australia
 PO Box 335, Pendle Hill, NSW 2145, Australia
 info@carrgroup.com.au

CAUTION!

To avoid premature degradation of the surface coating never use an abrasive cloth, sponge or other similar device on the surface. Carefully follow manufacturer instructions for recommended exposure time of cleaner on surface. Wipe off cleaning product residue with a clean, soft cloth after cleaning and allow the surface to dry completely before use.

IMPORTANT INFORMATION:

Cleaners/disinfectants can cause discolouration/bleaching on some colours. To avoid discolouration, always check on a small area first.





Information contained here is purely related to the compatibility of cleaners with **Carrflex**[®] materials and in no way relates to the efficacy of any cleaner for infection control or the safety of any cleaner for skin contact.

Testing info is specific to the products or formulations listed and in no way constitutes approval for other products from any brand listed without being tested.

All approvals are based on manufacturer recommended dilution with water unless otherwise specified and approval is void for any improper dilution.

If a particular cleaner/disinfectant is not listed, we advise to test the specific product in advance or to contact us to conduct testing.

CARR NEW ZEALAND LIMITED

 Tel: 64 9 274 0869
 6 Ra Ora Drive, East Tamaki, Auckland 2013, New Zealand
 PO Box 204110, Highbrook, Auckland 2161, New Zealand
 info@carrgroup.co.nz

CARR AUSTRALIA PTY LIMITED

 Tel: 61 2 9636 9525
 1 Spireton Place, Pendle Hill, NSW 2145, Australia
 PO Box 335, Pendle Hill, NSW 2145, Australia
 info@carrgroup.com.au