

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|--------------------------------|---|---|--------------------------------------|-------------|
| Ecolab® | A-456-II | | Minimum 64 to 1 dilution with water | 02-Dec-14 |
| Ecolab® | Actichlor Plus | Sodium Dichloroisochyanurate 30-60% | Minimum 100 to 1 dilution with water | 20-Nov-12 |
| Quantum Pacific Limited | Actimax | | According to Manufacturer Spec | 11-Apr-06 |
| ABC Compounding Co., Inc. | Aero Stain Away | | According to Manufacturer Spec | 11-Aug-16 |
| SCIENCECompany | Ammonium Chloride | Cocoalkyldimethylbenzyl ammonium chloride | According to Manufacturer Spec | 17-May-19 |
| Milton® | Anti-Bacterial Tablets | Sodium Dichloroisochyanurate 30-60% | 1 Tablet:2L Water | 19-Jul-12 |
| Ecolab® | Asepti Active Liquid | | According to Manufacturer Spec | 20-Nov-12 |
| Ecolab® | Asepti Active Wipes | Dual-Quaternary/Alcohol | Ready to Use | 20-Nov-12 |
| Synergy Health | Azomax™ Multi Surface Wipes | Benzylkoniumchloride 21.6% | Ready to Use | 15-Mar-13 |
| Synergy Health | Azowipe™ Hard Surface Wipes | Benzylkoniumchloride 2% | Ready to Use | 15-Mar-13 |
| PreVasive® | BAC Botanical Antimicrobial Cleaner | Thyme Oil | According to Manufacturer Spec | 22-May-19 |
| MEDCHEM® | BAK 1:750 | Benzalkonium Chloride | According to Manufacturer Spec | 12-May-19 |
| Clorox® | Bleach Wipes (Sodium Hypochlorite 0.1 - 1%) | Sodium Hypochlorite | According to Manufacturer Spec | 11-Aug-16 |
| Spartan Chemical Company Inc. | BNC-15™ | | According to Manufacturer Spec | 11-Aug-16 |
| Ameriderm Laboratories Limited | Care Source Disinfectant | | According to Manufacturer Spec | 10-Nov-11 |
| Metrex® Research | CaviCide® | Isopropanol | According to Manufacturer Spec | 06-Oct-16 |
| Metrex® Research | CaviCide1™ | Isopropanol | According to Manufacturer Spec | 05-Jan-17 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|---------------------------------|---|---|--------------------------------------|-------------|
| Metrex® Research | CaviWipes1™ | Isopropanol | Ready to Use | 05-Jan-17 |
| Metrex® Research | CaviWipes™ | Isopropanol | Ready to Use | 06-Oct-16 |
| Helix Solutions | Chlor-Clean | Sodium Dichloroisocyanurate 20-30% | 1000PPM (1 Tablet per 1 Litre Water) | 10-Oct-11 |
| Agar Cleaning Systems Pty. Ltd. | Chloradet | Sodium Hypochlorite <10% | Minimum 160 to 1 dilution with water | 06-Jun-16 |
| Ecolab® | Chlorwhite | Sodium Hypochlorite <10%, Sodium Hydroxide <10% | Minimum 100 to 1 dilution with water | 22-Oct-13 |
| Apollo Corporation | CID A-L™ II Quaternary Disinfectant | Quaternary Ammonium Compounds | According to Manufacturer Spec | 26-Sep-12 |
| Wexford Labs, Inc. | CleanCide® Wipes | Citric Acid | Ready to Use | 15-May-19 |
| Clorox Disinfecting Wipes | Clorox Disinfecting Wipes | | Ready to Use | 14-Feb-11 |
| Alco Mandala Products | DC-7 | | According to Manufacturer Spec | 16-Feb-12 |
| Delaval | Delaval Acid Wash | | Minimum 20 to 1 dilution with water | 06-Jun-14 |
| Clinell® | Detergent Wipes | Propane 1,2-diol | According to Manufacturer Spec | 27-May-09 |
| ArjoHuntleigh | Disinfectant Cleanser IV | | According to Manufacturer Spec | 23-Nov-09 |
| Clorox® | Dispatch® Disinfectant | | According to Manufacturer Spec | 16-Feb-12 |
| Tristel Solutions LTD | Distel Surface Disinfectant (Formerly Trigene®) | Polymeric Biguanide Hydrochloride <10% | According to Manufacturer Spec | 24-Oct-06 |
| Sealed Air Diversey Care | Divercleanse | Sodium Hypochlorite <5%, Sodium Peroxide <2% | Minimum 33 to 1 dilution with water | 15-Dec-09 |
| Buckeye® International, Inc. | Eco E23 S23 | | Minimum 256 to 1 dilution with water | 05-Jan-17 |
| Velocity Chemicals | Eco-Pine | | Minimum 20 to 1 dilution with water | 11-Apr-06 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|-------------------------------|--|---|-------------------------------------|-------------|
| Oxiver | Five 16 | | Minimum 16 to 1 dilution with water | 08-Nov-19 |
| Clorox HEALTHCARE® | Fuzion™ | | According to Manufacturer Spec | 05-Jan-17 |
| Spartan Chemical Company Inc. | GS Neutral Disinfectant #3502 | | According to Manufacturer Spec | 02-Dec-14 |
| Clinell® | Hand & Surface Sanitizer | Benzalkonium Chloride | According to Manufacturer Spec | 27-May-09 |
| 3M | HB Quat Disinfectant | Alkyl dimethyl benzyl ammonium chloride | Ready to Use | 06-Sep-19 |
| Spartan Chemical Company Inc. | hdq®C2 | Ammonium Chloride | According to Manufacturer Spec | 30-Jun-17 |
| Brighton Professional | Hepastat 256 | | According to Manufacturer Spec | 27-Aug-14 |
| Wipesplus | Hospital/Surgical Wipes | Ammonium Chloride | Ready to Use | 15-May-19 |
| (Various) | Hydrogen Peroxide 5% | | Ready to Use | 25-May-12 |
| Clorox HEALTHCARE® | Hydrogen Peroxide Cleaner Disinfectant Wipes | Hydrogen Peroxide | Ready to Use | 12-Apr-19 |
| Apollo Corporation | Hygena™ | Ethylalcohol 1.4% | According to Manufacturer Spec | 25-Sep-12 |
| Ecolab® | Incidin® Extra | Benzylchloride-Glucoprotamine 2% | According to Manufacturer Spec | 06-Jun-14 |
| Ecolab® | Incidin® Foam | Glucoprotamine 26% | According to Manufacturer Spec | 06-Jun-14 |
| Ecolab® | Incidin® Plus | Glucoprotamine 26% | According to Manufacturer Spec | 06-Jun-14 |
| (Various) | Isopropyl Alcohol 70% | | Ready to Use | 25-Feb-10 |
| Kimberly-Clark® | Isowipes® | Isopropanol 30-60% | Ready to Use | 19-Jul-12 |
| Schulke & Mayr GMBH | Kodan® Tinkture Forte | Hydroperoxide Solution 30% | According to Manufacturer Spec | 06-Jun-14 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|-------------------------------------|--|---|-------------------------------------|-------------|
| Legge Systems | Legcide HB Plus | | According to Manufacturer Spec | ***** |
| Whiteley Medical | Lemex | | According to Manufacturer Spec | 15-Mar-11 |
| Livingstone International (PTY) LTD | LIV-Wipes Jumbo | 70% Isopropyl Alcohol | Ready to Use | 15-Mar-11 |
| Sulco LTD | Mediwipes | | Ready to Use | 27-May-09 |
| Medline Industries Inc. | MICRO-KILL™ + Disinfecting, Deodorizing, Cleaning Wipes | Isopropyl Alcohol | According to Manufacturer Spec | 09-Mar-20 |
| Medline Industries Inc. | MICRO-KILL™ AF2 Wipes | Ammonium Chloride | Ready to Use | 15-May-19 |
| Medline Industries Inc. | MICRO-KILL™ Bleach Germicidal Wipes | Sodium Hypochlorite | Ready to Use | 15-May-19 |
| Medline Industries Inc. | MICRO-KILL™ Bleach-Free, Alcohol-Free Wipes | Ammonium Chloride | Ready to Use | 15-May-19 |
| Medline Industries Inc. | MICRO-KILL™ One Wipes | Ethyl Alcohol | Ready to Use | 15-May-19 |
| Medline Industries Inc. | MICRO-KILL™ Q3 Concentrated Disinfectant, Cleaner & Deodorizer | | According to Manufacturer Spec | 01-Feb-18 |
| Pro Restore® Products | Microban X-590 Institutional Spray | Diisobutylphenoxyethoxy ethyl dimethyl benzyl ammonium chloride | According to Manufacturer Spec | 11-Aug-16 |
| Anti-Slip Anywhere | Multi Kleen | | According to Manufacturer Spec | 15-Mar-11 |
| Medizar® | Multi-Surface Sanitiser Wipes | | Ready to Use | 17-May-19 |
| Ecolab® Cleantec | Nature Clean | | Minimum 38 to 1 dilution with water | 17-Dec-10 |
| Certo® | Non Acid Bathroom Cleaner | | According to Manufacturer Spec | 07-May-12 |
| PreVasive® | Noroxycdiff™ | Hydrogen Peroxide | According to Manufacturer Spec | 22-May-19 |
| Ecolab® | Oasis Pro 12 Neutral Cleaner | Alcohols <10% | Minimum 20 to 1 dilution with water | 06-Mar-12 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|--------------------------------|---|---|-------------------------------------|-------------|
| Schülke & Mayr GmbH | Ocetenisept® | 100g solution contains: Octenidine dihydrochloride 0.1g, 2-Phenoxyethanol | According to Manufacturer Spec | 06-Jun-14 |
| Whiteley Medical | Odacon Spray | | According to Manufacturer Spec | 14-Jun-11 |
| OdorXit Naturals | Odor Eliminator Spray | | Ready to Use | 22-Feb-18 |
| Oxiver | One Wipes | | Ready to Use | 14-Feb-20 |
| Sealed Air Diversey Care | Oxivir® Five | | Minimum 16 to 1 dilution with water | 08-Nov-19 |
| Sealed Air Diversey Care | Oxivir® TB Liquid and Wipes | | Ready to Use | 22-Apr-13 |
| Ecolab® | Oxycide | | According to Manufacturer Spec | 06-Jul-14 |
| Pal® International | PAL® Detergent Wipes | | Ready to Use | 27-May-09 |
| Pal® International | PAL® Disinfectant Wipes | | Ready to Use | 27-May-09 |
| PCS Process Cleaning Solutions | PCS 7000 | Sodium Hypochlorite | According to Manufacturer Spec | 06-Oct-17 |
| PCS Process Cleaning Solutions | PCS Spray & Wipe Hand and Surface Cleaner | | Ready to Use | 06-Oct-17 |
| Pure Green, LLC | Puregreen 24 (Pure Green) | | According to Manufacturer Spec | 30-Jan-09 |
| Reynard Health Supplies | Reynard Neutral Detergent Wipes | | Ready to Use | 08-Feb-07 |
| Reynard Health Supplies | Reynard Surface Sanitizing Wipes | | Ready to Use | 27-May-11 |
| Sealed Air Diversey Care | Room Care R2L | | According to Manufacturer Spec | 24-Oct-12 |
| STERI-7® | S-7 Xtra Wipes | | Ready to Use | 28-Nov-16 |
| Buckeye® International, Inc. | Sanicare TBX™ | | According to Manufacturer Spec | 21-Sep-10 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|-------------------------------|--|------------------------------------|--------------------------------------|-------------|
| Buckeye® International, Inc. | Sanicare™ Disinfectant Wipes | | Ready to Use | 21-Sep-10 |
| Alco Mandala Products | Sanni-Rinse | | According to Manufacturer Spec | 16-Feb-12 |
| (Various) | Sodium Hypochlorite | | Minimum 77 to 1 dilution with water | ***** |
| Hu-Friedy Mfg. Co., LLC. | Sporeclear RTU Spray or Wipes | Benzalkonium Chloride 0-0.5% | Minimum 13 to 1 dilution with water | 31-Jan-16 |
| Clinell® | Sporicidal Wipes | | Ready to Use | 27-May-09 |
| Ecolab® | Spring Clean | | Minimum 20 to 1 dilution with water | 11-Apr-06 |
| Noble Pine Products Co. | Steri-Fab® | Isopropyl Alcohol | According to Manufacturer Spec | 11-Aug-16 |
| Hartmann | Sterillium® | Glucoprotamine 30% | According to Manufacturer Spec | 06-Jun-14 |
| sBioMed, LLC | Steriplex SD Disinfectant Sporicide Wipes (Without Bleach) | Hydrogen Peroxide | According to Manufacturer Spec | 20-May-15 |
| Sterri-Matt® | Sterri-Matt® (Pure) Foam | | Ready to Use | 15-Dec-09 |
| Sealed Air Diversey Care | Stride® Citrus Smartdose | Ethoxylated Alcohol 20-30% | Minimum 750 to 1 dilution with water | 01-Jun-15 |
| Sealed Air Diversey Care | Stride® HC (JF) | Ethoxylated Alcohol 20-30% | According to Manufacturer Spec | 01-Jun-15 |
| Sealed Air Diversey Care | Suma San Conc D4A | | According to Manufacturer Spec | 28-Oct-14 |
| Sealed Air Diversey Care | Suma San Conc D4A | Didecyl Dimethyl Ammonium Chloride | According to Manufacturer Spec | 28-Oct-14 |
| Spartan Chemical Company Inc. | Super HDQ Neutral® | | Minimum 256 to 1 dilution with water | 11-Aug-16 |
| PDI | Super Sani Cloth® Germicidal Disposable Wipes | | Ready to Use | 02-Sep-10 |
| Anios Laboratories | Surfa 'Safe | Chlorhydrate 1.4% | According to Manufacturer Spec | 06-Jun-14 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|------------------------------|-----------------------------------|---|--------------------------------------|-------------|
| Bioscience Fresche® | Surface Plus | | According to Manufacturer Spec | 09-Feb-15 |
| Reynard Health Supplies | Surface Sanitising Wipes | | Ready to Use | 18-Feb-09 |
| 3M | TB Quat Disinfectant | Ammonium Chloride | Ready to Use | 09-Mar-18 |
| Buckeye® International, Inc. | Terminator™ One Step Disinfectant | | According to Manufacturer Spec | 21-Sep-10 |
| Apollo Corporation | Therasol™ | | According to Manufacturer Spec | 25-Sep-12 |
| Medline Industries Inc. | Touch Screen Cleaning Wipes | 70% Isopropyl Alcohol | Ready to Use | 15-May-19 |
| Hunters Products | TRI-DET | | Minimum 80 to 1 dilution with water | 02-Dec-14 |
| (Various) | Trix (Soap) | | Ready to Use | 27-Apr-10 |
| Verna Care | Tuffie 5 Cleaning Wipes | Cocoalkyl Dimthybenzyl Ammonium Chloride 0-0.5% | Ready to Use | 19-Jul-12 |
| Apollo Corporation | Turbo Clean™ | Ethanol 9.5% | According to Manufacturer Spec | 25-Sep-12 |
| Texwipe® | TX®3270 | 70% isopropyl Alcohol | Ready to Use | 12-Apr-19 |
| Clinell® | Universal Disinfectant Wipes | | According to Manufacturer Spec | 27-May-09 |
| Clinell® | Universal Wipes | | Ready to Use | 12-May-19 |
| True Blue Chemicals | Use All | Glycol 1-6% | According to Manufacturer Spec | 28-Sep-12 |
| Whiteley Medical | Viraclean® | Benzalkoniumchloride | According to Manufacturer Spec | 19-Jul-12 |
| Sealed Air Diversey Care | Virex® II 256 | | Minimum 256 to 1 dilution with water | 07-May-12 |
| Sealed Air Diversey Care | Virex® Plus | Ammonium Chloride | According to Manufacturer Spec | 19-Dec-18 |

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

Sorted by Product Name:

| MANUFACTURER | PRODUCT | MAIN ACTIVE COMPONENT | CONCENTRATION | DATE TESTED |
|--------------|-------------|-----------------------|--------------------------------------|-------------|
| Dupont™ | Virkon® | | Minimum 20 to 1 dilution with water | 15-Dec-09 |
| Wexford Labs | Wexcide 128 | | Minimum 128 to 1 dilution with water | 11-Aug-16 |

CAUTION!

To avoid premature degradation of the surface coating never use an abrasive cloth, sponge or other similar device on the surface. Carefully follow manufacturer instructions for recommended exposure time of cleaner on surface. Wipe off cleaning product residue with a clean, soft cloth after cleaning and allow the surface to dry completely before use.

IMPORTANT INFORMATION:

Cleaners/disinfectants can cause discolouration/bleaching on some colours. To avoid discolouration, always check on a small area first.

Information contained here is purely related to the compatibility of cleaners with Carrflex® materials and in no way relates to the efficacy of any cleaner for infection control or the safety of any cleaner for skin contact.

Testing info is specific to the products or formulations listed and in no way constitutes approval for other products from any brand listed without being tested.

All approvals are based on manufacturer recommended dilution with water unless otherwise specified and approval is void for any improper dilution.

If a particular cleaner/disinfectant is not listed, we advise to test the specific product in advance or to contact us to conduct testing.

Listed chlorine based cleaners should be used only in instances of irregular spot cleaning and bodily fluid spills.

CARR NEW ZEALAND LIMITED

Tel: 64 9 274 0869
6 Ra Ora Drive, East Tamaki, Auckland 2013, New Zealand
PO Box 204110, Highbrook, Auckland 2161, New Zealand
info@carrgroup.co.nz

CARR AUSTRALIA PTY LIMITED

Tel: 61 2 9636 9525
1 Spireton Place, Pendle Hill, NSW 2145, Australia
PO Box 335, Pendle Hill, NSW 2145, Australia
info@carrgroup.com.au